Intended Use: The CONTOUR[®]NEXT test strips are intended for use by persons with diabetes for self-testing, and by health care professionals for use on a single patient. The CONTOUR NEXT test strips are intended to be used to monitor glucose concentrations in whole blood.

Summary: The CONTOUR NEXT test strip is designed for use with all available CONTOUR[®]NEXT blood glucose meters. The test provides a quantitative measurement of glucose in blood for CONTOUR NEXT meters. The available CONTOUR NEXT systems are not intended for diagnosis of, or screening for diabetes mellitus.

Storage and Handling:

- Store the strips at temperatures between 0°C and 30°C.
- **Store test strips in their original bottle only. Always close the lid immediately and tightly after removing a test strip.**
- **Wash and dry your hands well before handling the test strips.**
- Do not use the test strips after the expiry date. The expiry date is printed on the bottle label and on the test strip carton.
- If the meter and/or test strips are moved from one temperature to another, allow 20 minutes for them to adjust to the new temperature before performing a blood glucose test. Your user guide will identify the appropriate operating temperature range for the CONTOUR NEXT meter you are using.
- The test strips are for single use only. **Do not reuse test strips.**
- Upon opening the test strip carton ensure the strip bottle cap is securely closed. If the cap is not closed, do not use the test strips for testing. Examine the product for missing, damaged, or broken parts. Contact Customer Service at 1-800-268-7200 for replacement parts and assistance.

Test Procedure: See your CONTOUR NEXT user guide and accompanying inserts for detailed testing instructions.

Test Results:

Your meter has been preset to display results in mmol/L (millimoles of glucose per liter). Results in mmol/L will **always** have a decimal point (e.g., 5.3 $\frac{\text{mmol}}{\text{L}}$); results in mg/dL will **never** have a decimal point (e.g., 96 $\frac{\text{mg}}{\text{dL}}$). If your test result is displayed in mg/dL, contact Customer Service at 1-800-268-7200.

- If your blood glucose reading is under 2.8 mmol/L on the meter display, follow medical advice immediately.
- If your blood glucose reading is above 13.9 mmol/L on the meter display, call your health care professional as soon as possible.
- Always consult your health care professional before changing your medication based on CONTOUR NEXT test results.

Blood glucose values will vary from time to time depending on food intake, medication dosages, health, stress or exercise. Consult your health care professional for the target value appropriate for you. Standard medical practice goals for a typical non-pregnant individual with diabetes are¹

- Fasting Glucose: 4.0 to 7.0 mmol/L
- Postprandial Glucose: (1-2 hours after meal) 5.0 to 10.0 mmol/L

Questionable or Inconsistent Results: See the CONTOUR NEXT user guide for problem solving. If attempts to correct a problem fail, contact Customer Service at 1-800-268-7200.

Quality Control: You should perform a control test when using your meter for the first time, or when you open a new bottle or package of test strips, or if you think your meter may not be working properly, or if you have repeated unexpected blood glucose results. Only use CONTOUR[®]NEXT control solutions. These control solutions are designed specifically for use with the CONTOUR NEXT systems. The control results should fall within the control range(s) printed on the test strip bottle label and the test strip carton. If they don't, do not use your meter for blood glucose testing until you resolve the issue.

CAUTION

Information For Safety

- For *in vitro* diagnostic use only. External use, do not swallow.
- Potential Biohazard: All products or objects that come in contact with human blood, even after cleaning, should be handled as if capable of transmitting infectious diseases. The user should follow the recommendations for the prevention of blood-borne transmissible diseases in health care settings as recommended for potentially infectious human specimens.²
- Dispose of the test strips as medical waste or as advised by your health care professional.

Chemical Composition: FAD glucose dehydrogenase (*Aspergillus sp.*, 4.0 U/test strip) 21%; phenothiazine mediator 54%; nonreactive ingredients 25%.

Principles of the Procedure: The CONTOUR NEXT blood glucose test is based on measurement of electrical current caused by the reaction of glucose in the blood with the reagents on the electrode of the strip. The blood sample is drawn into the tip of the test strip through capillary action. Glucose in the sample reacts with FAD glucose dehydrogenase (FAD-GDH) and a phenothiazine derivative as a mediator. Electrons are generated, producing a current that is proportional to the glucose in the sample. After the reaction time, the glucose concentration in the sample is displayed. No calculation is required by the user.

Comparison Options: Available CONTOUR NEXT systems are designed for use with venous and capillary whole blood. Comparison to a laboratory method must be done simultaneously with aliquots of the same sample. Note: Glucose concentrations drop rapidly due to glycolysis (approximately 5%–7% per hour).³

Limitations:

1. **Preservatives:** Blood may be collected by health care professionals into test tubes containing heparin. Do not use other anticoagulants or preservatives.
2. **Altitude:** Up to 6301 meters does not significantly affect results.
3. **Alternative Site Testing:** Please see your meter user guide for Alternative Site Testing instructions.
4. **Peritoneal dialysis solutions:** Icodextrin does not interfere with CONTOUR NEXT test strips.
5. **Contraindications:** Capillary blood glucose testing may not be clinically appropriate for persons with reduced peripheral blood flow. Shock, severe hypotension, hyperosmolar hyperglycemia and severe dehydration are examples of clinical conditions that may adversely affect the measurement of glucose in peripheral blood.⁴
6. **Interference:** The CONTOUR NEXT system was tested against the following potentially interfering substances occurring naturally in the blood: bilirubin, cholesterol, creatinine, galactose, glutathione, hemoglobin, triglycerides, and uric acid. No interfering effect was observed for any substance at the highest concentration⁵ of either the common pathological level or three times the upper reference value.⁶
7. **Interference:** The CONTOUR NEXT system was tested against the following potentially interfering substances occurring from therapeutic treatments: ascorbic acid, paracetamol (acetaminophen), dopamine, sodium gentisate, ibuprofen, icodextrin, L-dopa, maltose, methyl dopa, pralidoxime iodide, sodium salicylate, tolazamide, tolbutamide. No interfering effect was observed for any substance at the highest concentration⁵ of either the toxic concentration or three times the maximum therapeutic concentration.⁶
8. **Xylose:** Do not use during or soon after xylose absorption testing. Xylose in the blood will cause interference.
9. **Hematocrit:** CONTOUR NEXT test strip results are not significantly affected by hematocrit levels in the range of 0% to 70%.⁶

See your meter user guide for intended use with neonates and alternative site testing.

Technical Information:

Please refer to your CONTOUR NEXT meter user guide for accuracy and precision information.

Performance Evaluations (for health care professionals): Please contact Customer Service at 1-800-268-7200 for data and results related to the performance of your CONTOUR NEXT system.

References:

1. Canadian Diabetes Association 2013 Clinical Practice Guidelines for the Prevention and Management of Diabetes in Canada. *Canadian Journal of Diabetes*. 2013;37(supplement 1):S33.
2. Sewell DL. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline, 3rd Edition*. Clinical and Laboratory Standards Institute. CLSI document M29-A3; ISBN 156238-567-4. March 2005.
3. Burtis CA, Ashwood ER, editors. *Tietz Fundamentals of Clinical Chemistry*. 5th edition. Philadelphia, PA: WB Saunders Co; 2001;444.
4. Atkin SH, et al. Fingertick glucose determination in shock. *Annals of Internal Medicine*. 1991;114(12):1020-1024.
5. McEnroe RJ, et al. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. EP7-A2, vol 25, no 27. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
6. Bernstein R, Parkes JL, Goldy A, et al. A new test strip technology platform for self-monitoring of blood glucose. *Journal of Diabetes Science and Technology*. 2013;7(5):1386-1399.

Manufactured for:
Ascensia Diabetes Care Holdings AG
Peter Merian-Strasse 90
4052 Basel, Switzerland
1-800-268-7200
www.diabetes.ascensia.com

For information on patents, and related licenses, see:
www.patents.ascensia.com

Ascensia, the Ascensia Diabetes Care logo, Clinilog, Contour, Glucofacts, and the No Coding logo are trademarks and/or registered trademarks of Ascensia Diabetes Care Holdings AG.

TM see www.ascensidiabetes.ca/en/tm-mc

PROCESS K		FOLD	
RR DONNELLEY 5050 1st Ave South Suite #101 Seattle, WA 98134 Tel: 206.587.0800 Fax: 206.587.0356	FILE NAME: 85306324_CntrNXT_Rgtns_EN_FRca.indd		NEW COMP BAN: 85306324
	DESCRIPTION: Contour Next Reagent Insert		REV: 05/16
	DATE: 07/13/16	CLIENT: Ascensia Diabetes Care	V0 06/16/16 Initial Build
	PROGRAM VERS: InDesign CC(2014)	SKU: Multi	V1 06/24/16 DMG Annotations
			V2 07/07/16 DMG Annotations
		V3 07/13/16 DMG Annotations	
MARKET: CA	SPEC: 1000008285 REV 05	BAR CODE: 85306324 (Code 128 & 2D Matrix)	
LANGUAGES: EN/FRca	PARENT BAN: Multi		
TRIM SIZE: 11.654" (H) x 6.614" (W)	PREVIOUS BAN: 84603759		
	MASTER REF: XXXXX		
LASER PROOF	EPSONS PROOF	SIGNATURE	VERSION: V3
<input type="radio"/> APPROVED	<input type="radio"/> APPROVED	PROJ. MGR. SIGNATURE _____	DATE _____
<input type="radio"/> RESUBMIT WITH CHANGES	<input type="radio"/> RESUBMIT WITH CHANGES	PROJ. MGR. SIGNATURE _____	DATE _____
		PROJ. MGR. SIGNATURE _____	DATE _____

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Références :

1. Canadian Diabetes Association 2013 Clinical Practice Guidelines for the Prevention and Management of Diabetes in Canada. *Canadian Journal of Diabetes*. 2013;37(supplément 1):S33.
2. Sewell DL. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline, 3rd Edition*. Clinical and Laboratory Standards Institute. CLSI document M29-A3; ISBN 156238-567-4. March 2005.
3. Burtis CA, Ashwood ER, editors. *Tietz Fundamentals of Clinical Chemistry*. 5th edition. Philadelphia, PA: WB Saunders Co; 2001;444.
4. Atkin SH, et al. Fingerstick glucose determination in shock. *Annals of Internal Medicine*. 1991;114(12):1020-1024.
5. McEnroe RJ, et al. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. EP7-A2, vol 25, no 27. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
6. Bernstein R, Parkes JL, Goldy A, et al. A new test strip technology platform for self-monitoring of blood glucose. *Journal of Diabetes Science and Technology*. 2013;7(5):1386-1399.

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